

# Hand Engine

---



35,000 RPM motor

Autoclavable hand piece and cord

Rotation easily switched from forward to reverse.

Self locking head - No more twisting and turning to remove. Just pull apart.

Use either manual speed control knob or variable speed foot switch.

Many accessories available

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Diamond Fraises

## Diamond Fraises - for dermabrasion & cutaneous prophylaxis

The largest selection of sizes and grades of abrasion available anywhere

Actual Sizes



	Regular	Coarse	Extra Coarse
17x14mm.	\$65	\$124	\$145
17x10mm.	52	110	131
17x8mm.	44	94	111
15x10mm.	52	94	111
17x6mm.	37	79	93
17x4mm.	29	64	74
17x2.5mm.	24	47	57
13x2.5mm.	21	44	
10x2.5mm.	21	42	
CONE	23	30	
PEAR	23	30	37
DOME	91	123	160
BULLET	40	53	69
FRUSTUM	68	108	140
FLAME	32	43	52
BUTTON	34	46	59
MINI BULLET	26	36	53



**Top Row:** 10x2.5mm, Cone

**Bottom Row:** Button, Mini-Bullet

- one piece construction for 100% centricity


- diamond chips grade for grade are the largest and provide for maximum abrasion

- all edges have been rounded, avoiding dangerous sharp corners

## Stainless Steel Wire Planing Brushes

All fraises and brushes have 3/32" diameter shaft which is the international standard

Actual Sizes



Small - 17x3mm. (regular, soft, course) \$69

Medium - 17x6mm. 72

Large - 17x9mm. 79

Extra Large - 17x12mm. 82

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)