

Liposuction Aspirator



Our liposuction aspirator features:

- Advanced state of the art piston pumps
- Maximum deep vacuum quickly reached and maintained
- Dual headed pumps eliminating back pressure
- Pumps set in series rather than in tandem
- Three-quarter horsepower
- Quiet, vibration-free, reliable operation
- Standard implosion-proof canisters and disposable liners
- Five-year warranty on all parts and labor

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

cosmeticsurgerysuppliers.com

Aspirator Accessories



Standard Suction Tubing



Flexi-Tubing



Disposable HEPA Filter



2000cc Abbott Canister
with Disposable Liner



1500cc Baxter Canister
with Disposable Liner

Purchase from:





Cosmetic Surgery Suppliers, Inc.
7758 Hampton Place
Loganville, GA 30052
800.525.7752
770.934.0235

cosmeticsurgerysuppliers.com


Syringe Aspiration

syringe cannula for toomey tip
(made to order)



toomey tip syringe

syringe-to-syringe transfers



Syringe Locks



snapper type

- 10 cc
- 30 cc
- 60 cc



horseshoe type

- 10 cc
- 30 cc
- 60 cc

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

cosmeticsurgerysuppliers.com

Disposable Tubing



Standard Tubing

- 3/8" I.D.
- Transparent for easy visualization
- 10 foot lengths for easy use
- 10 pieces to a case



Flexi-Tubing

- 9/32" I.D.
- Super flexible tubing
- 12 foot lengths for easy use
- 10 pieces to a case

Tumescent Infiltration Tubing



AFT Tubing

- Bag spike and luer-lock connector
- Thicker section is silicone, to avoid latex reactions
- 10 pieces to a case



T.I. Tubing

- Bag spike and luer-lock connector
- 8 foot lengths for convenience
- 10 pieces to a case

Purchase from:



Cosmetic Surgery Suppliers, Inc.
7758 Hampton Place
Loganville, GA 30052
800.525.7752
770.934.0235

cosmeticsurgerysuppliers.com

Cannula & Brushes

Cosmetic Surgery Suppliers, Inc. Cannulas are custom made to your specifications when you place your order. When ordering your custom cannula please specify the following four items: tip style, lumen diameter, lumen length and handle style. All of the most popular tip styles are available and custom tip styles can also be manufactured.



One-Hole Standard



Candy Cane



Three-Hole Mercedes



Spatula




Three-Hole Accelerator




Three-Hole Cage

Standard Diameters
2mm, 3mm, 4mm, 5mm & 6mm

Standard Lengths
10cm, 15cm, 23cm & 30cm



One-piece cannula handle



Luer lock hub



60cc syringe hub

Cannula Cleaning Brushes

Available in 5 sizes (2mm, 3mm, 4mm, 5mm & 6mm), the brushes fit the full range of cannula sizes currently being used in most surgeries



Cosmetic Surgery Suppliers Inc. Cannula and tumescent needles are manufactured to the highest industry standards for quality. They include a two year warranty on materials & workmanship

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

cosmeticsurgerysuppliers.com